

## SUPPLEMENTARY TABLE

Supplementary Table 1. Primers used for RQ-PCR, RQ-MSP and BSP.

	Primer sequences
RQ-PCR	
MiR-335-Forward	5'-TTATAAAGCAATGAGACTGATT-3'
MiR-335-Reverse	Manufacturer-provided miScript universal primer
U6-Forward	5'-GTGCTCGCTTCGGCAGCACATATAC-3'
U6-Reverse	5'-AAAATATGGAACGCTTCACGAATTTG-3'
ID4-Forward	5'-CATCCCGCCCAACAAGAAAGTCA-3'
ID4-Reverse	5'-GCCGGGTCGGTGTGAGCGCAGT-3'
ABL-Forward	5'-TCCTCCAGCTGTTATCTGGAAGA-3'
ABL-Reverse	5'-TCCAACGAGCGGCTTCAC-3'
RQ-MSP	
MiR-335-Methylation-Forward	5'-TCGTGTTTTTGGTGGTTATC-3'
MiR-335-Methylation-Reverse	5'-TAACCGCTCGAACTAAAAAAC-3'
ALU-Forward	5'-TTAGGTATAGTGGTTTATATTTGTAATTTTAGTA-3'
ALU-Reverse	5'-ATTAATAAATAATCTTAAACTCCTAACCTCA-3'
BSP	
MiR-335-Bisulfited-Forward	5'-GGGGTTTTGTTTTTTAATTGT-3'
MiR-335-Bisulfited-Reverse	5'-AAATACCCCAACTCTTCCTTA-3'

RQ-PCR: real-time quantitative PCR; RQ-MSP: real-time quantitative methylation-specific PCR; BSP: bisulfite sequencing PCR.