

SUPPLEMENTARY TABLES

Supplementary Table 1. Information of antibodies of CD8, B7-H3, and CD47.

Antibody	Clone	Source	Antigen retrieval condition	Antibody dilution	Antibody incubation time
CD8	C8/144B	Cell Signaling Technology	EDTA buffer, pH9.0	1:200	18hrs
B7-H3 (CD276)	D9M2L	Cell Signaling Technology	EDTA buffer, pH9.0	1:300	18hrs
CD47	D307P	Cell Signaling Technology	EDTA buffer, pH 9.0	1:300	18hrs

Please browse Full Text version to see the data of Supplementary Table 2.

Supplementary Table 2. The correlation coefficients between Th1/IFN γ signature genes and B7-H3 or CD47.