

SUPPLEMENTARY TABLE

Supplementary Table 1. Primers used for quantitative real time PCR.

Primer	Sequence
CRART16 Forward	5'-TGATAGTGAGGCCTCCTGCAA-3'
CRART16 Reverse	5'-CTGGAGTTCTGCAGGTTCCCTTT-3'
miR-122-5p Forward	5'-GCTGTGGAGTGTGACAATGGTGTGTTG-3'
U6 Forward	5'-GCAAGGATGACACGCAAATTC-3'
FOS Forward	5'-CAAGCGGAGACAGACCAACT-3'
FOS Reverse	5'-TGAGCCGCTAGGATGAACTC-3'
VEGFD Forward	5'-ATCTGTATGAACACCAGCACCTC-3'
VEGFD Reverse	5'-TGGCAACTTTAACAGGCACTAAT-3'
GAPDH Forward	5'-ATGGGGAAGGTGAAGGTCGG-3'
GAPDH Reverse	5'-GACGGTGCCATGGAATTTGC-3'